

Visit us at AUA 2017 Exhibit Booth Space 2112

The Problem: Ball tip not withstanding laser lithotripsy

NEW BALL TIP

USED BALL TIP




Ball tip has intense tip degradation

The Solution: Scopesafe Fiber with Soft Tip™

NEW SOFT TIP

USED SOFT TIP





Soft Tip™ has NO significant tip degradation

- √ No significant tip degradation
- Multiple passes through scope
- √ No need to strip fiber
- √Pin stone safely
- √ Reduces mucosal trauma
- √ Reduces digital interference

The Soft Tip Succeeds Where the Ball Tip Fails.

ScopeSafe laser fiber with Soft Tip is a trademark of Optical Integrity. Patent pending.

19408 Panama City Beach Pkwy • Panama City Beach • FL • 32413
850-233-5512 • www.sales@opticalintegrity.com

