

**The NewSurg KTP Affordably Replaces  
Old Iridex DioLite / VariLite and  
Laserscope Aura Lasers.**

**If you are frustrated with your Pulsed  
Dye Laser (expensive dye packs,  
service calls, purpura ...) the KTP Laser  
is your solution.**

According to David J. Goldberg, MD, JD with offices in New York and New Jersey, "This is a wonderful system for treating facial telangiectases, erythema, red scars and other vascular disorders. The KTP is effective and very easy to use and can be effortlessly moved from room to room. It requires no warm-up time, no disposables and no expensive warranty. Our NewSurg KTP delivers what I need in a dermatology laser system."



Richard Goldfarb, MD in Pennsylvania states, "The KTP's simplicity and ease-of-use make it a workhouse in our practice. I can precisely target a lesion with the desired energy dosage and effectively treat an older patient with pigmented lesions and sun damage and then quickly adjust to the exact settings to safely eliminate a pediatric hemangioma. With the KTP's Scan Mode, we treat large areas quickly and this is sometimes a better treatment choice than our intense pulsed light (IPL) system. We no longer need a pulsed-dye laser (PDL)."

For Ayman El Attar, MD, with multiple offices in New Jersey, the NewSurg KTP reliability targets melanin rich lesions on darker skin types, including dermatosis papulosa nigra (DPNs) with high specificity and virtually no hyper or hypopigmentation post treatment. "I treat a lot of ethnic skin and in contrast to electrocautery or cryo, which leaves patients discolored. I use the precision of my KTP to address pigmented lesions on darker skin types" explained El Attar. "The platform uses preset protocols for common lesions, and the lightweight pen-like handpiece provides predictable results with minimal collateral heating. Patients see immediate improvements to their lesions" continues Dr. El Attar.

"The KTP is also very reliable and we safely move our KTP from office to office in the rugged carry case which is included."



Owners of the NewSurg KTP appreciate the preset treatment parameters that make it safe and easy to use both novice and experienced physicians and technicians. It is the perfect dermatology laser!

To learn more about the most advanced and affordable dermatology laser for treating vascular and pigmented lesions, visit [www.newsurg.com](http://www.newsurg.com)

# NewSurg

Email [sales@newsurg.com](mailto:sales@newsurg.com) or call  
(215) 570-4327 with any questions  
or to schedule a demo